



**THE FUTURE OF LONG-TERM CARE:  
ECONOMICAL AND COMFORTABLE**

## **PB 325 III**

- ✓ Modern and compact design
- ✓ Low noise and reliable LINAK control
- ✓ 3-year warranty
- ✓ Easy transportation and assembly



Safe working load:  
**200 kg**

Height adjustment:  
**40 - 80 cm (standard)**  
**23 - 63 cm (option\*)**  
\*concerns low version

### New LINAK drive system

- ✓ compact control box CA20
- ✓ low noise electric actuators LA24,
- ✓ Handset with comfort function

**Lifting pole with triangle handle (standard)**

**Wooden front panel overlays (option)**



Side rails and new designed sliders allow to use mattress up to **15 cm**

# OPTIONS AND ADDITIONAL EQUIPMENT



**Metal full-length side rails /  
Metal slats**



**Upholstered cover for side  
rails**



**Bed lever S20**



**Bed extension 20 cm**



**Extra-height safety side rails**

# TECHNICAL SPECIFICATION OF BASIC CONFIGURATION

Height adjustment range

40 cm ÷ 80 cm

Safe working load

200 kg

External dimensions  
width x length x height

101 cm x 214 cm x 86 [126] cm

Mattress dimensions

90 cm x 200 cm x 15 cm

Mattress base

Wooden spring slats

Drive system control  
box/actuators/handsets

LINAK  
(CA20 / LA24 / HL74 Comfort)

Castors

100 mm (diameter), single  
break

Total weight Set of side rails  
/ Lifting pole

85,6 kg  
(11,2 kg / 5,3 kg)



Elbur sp. z o.o. sp.k.  
Działosza 34  
56-500 Syców,  
T: +48 62 786 97 80  
E: info@elbur.eu  
W: www.elbur.eu

Distributor