

PB 326



Our versatile PB 326 bed with its attractive design blends easily into many environments. Comfort for the user through a variety of functions, as well as ergonomic working conditions for the carer, make this product especially valued by rehabilitation equipment rental companies and individual users.



In the picture: PB 326 bed, standard design in beech. Accessories: cabinet Rubens 3, beech.

PB 326



Standard design

DETAILS:

- | | |
|---|----------------------------|
| 1. The beechwood elements | 4. Knee section adjustment |
| 2. Mattress support platform with firmness adjustment | 5. Brake castors |
| 3. Handset with lockable functions | 6. Lifting pole |


OPTIONS:

length and mattress support material choices • Trendelenburg position


BASIC TECHNICAL DATA:

Safe working load:	175,0 kg
Maximum user weight:	140,0 kg
Mattress support platform height adjustment	from 40 to 80 cm
Back rest adjustment:	0 ÷ 70°
Lower leg rest adjustment:	0 ÷ 20°
External dimensions (length × width):	216 × 104 cm
Clearance under the bed:	approx. 22 cm
Total weight:	90,0 kg
Castors:	100 mm with a brake


FUNCTIONS:



ADJUSTABLE LEG SECTIONS



ADJUSTABLE BACK SECTIONS



ADJUSTABLE HEIGHT OF MATTRESS SUPPORT PLATFORM



ANTI-TRENDELENBURG POSITION