

# PB 536




Modern design and multi-functionality are the dominant features of the PB 536 bed. Increased working load up to 200 kg, hidden castors and central brake facilitate the ease of work for the staff which are proven solutions used in many nursing homes and increasing more often by individual clients within a home care setting.



In the photo: The bed PB 536 beech, with divided rails on one side of the bed. Accessories: Rubens 9, beech.

# PB 536



PB 536 bed with divided railings on one side

## DETAILS:

- |   |                                    |
|---|------------------------------------|
| 1. Maximum safe working load 200 kg                           | 4. Handset with lockable functions |
| 2. Central brakes   | 5. Hidden castors                  |
| 3. Latest technology of lifting the mattress support platform | 6. Divided side rails (option)     |

## OPTIONS:

mattress support platform length choice • colour and mattress support material choices • divided side rails


## BASIC TECHNICAL DATA:

Safe working load:	200,0 kg
Maximum user weight:	165,0 kg
Mattress support platform height adjustment	from 39 to 80 cm
Back rest adjustment:	0 ÷ 70°
Lower leg rest adjustment:	0 ÷ 20°
External dimensions (length × width):	207 × 106 cm
Clearance under the bed:	approx. 16 cm
Total weight:	131,2 kg
Hidden castors:	50 mm


## FUNCTIONS:



ADJUSTABLE LEG SECTIONS



ADJUSTABLE BACK SECTIONS



ADJUSTABLE HEIGHT MATTRESS SUPPORT PLATFORM